

User manual  
**Telephone-Kit Comfort Contego**

Works with most telephones which have a modular connector to the handset

- Components:** 1) Switch box with telephone cable    2) Cable (1)    3) Cable (2)




**Getting started**


1. Switch on your Contego transmitter and receiver by pressing the On/Off-button for approx. 2 seconds.
2. Enter the menu of the Contego transmitter (T900). Choose *Input - Microphone*. A small microphone icon shows in the down left corner of the display.
3. Connect the induction loop to your Contego receiver (R900) and turn the hearing aid to the T-position. If you do not use hearing aids, you can use headphones instead.

**Connect the Telephone-Kit to the Comfort Contego**


1. Disconnect the handset cable at the telephone.



2. I) Connect the switch box cable to the telephone.  
 II) Connect the handset cable to the switch box.  
 III) Connect the cable (1) to the switch box and the line in of the Contego transmitter T900.



3. I) Pick up the handset and check if you can hear the dial tone (via your headphones or the induction loop).  
 II) If you cannot hear the dial tone, try cable(2); now you should be able to hear the dial tone.
4. The telephone volume can be adjusted with the arrows and at the Contego transmitter (T900).



**Make / Receive a call**

Use the telephone as normal. You talk into the handset.

The voice of your communication partner is sent from the Contego transmitter T900 to your Contego receiver R900.

You hear with headphones or via the induction loop with your hearing aid in T-position.

You can adjust the volume of the call with the volume control of your Contego receiver (R900) during an ongoing call.



**COMFORT  
CONTEGO**

A quality product from Sweden



**COMFORT  
CONTEGO**

A quality product from Sweden

## TELEPHONE-KIT

User manual

### Telephone-Kit Comfort Contego

Works with most telephones which have  
a modular connector to the handset



For additional information, please see the User Guide  
for Comfort Contego or our website [www.comfortaudio.us](http://www.comfortaudio.us)



**Comfort Audio Inc**

350 South Northwest Hwy, Ste 300

Park Ridge, IL 60068 USA

Phone (847) 656 5450

[usa@comfortaudio.com](mailto:usa@comfortaudio.com) [www.comfortaudio.us](http://www.comfortaudio.us)

MA0336